


# MANUAL

## **BIX-LV39 ADVANCED DIGITAL RECTAL EXAMINATION MODEL**



**Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd.**

**Ada Med Supply Limited**

## **BIX-LV39 ADVANCED DIGITAL RECTAL EXAMINATION MODEL**

### **Product Features:**

1. Digital rectal examination is the must and the simplest method for benign prostatic hyperplasia (BPH) and rectal cancer check.
2. Model anatomical structure clearly, it has the anus, rectum and prostate and other structures.
3. The model is vivid, feel real, it is easy for replacement connect components of prostate and rectum.
4. This product is made of imported plastic materials, stainless steel die cast molding processes, materials refinement, structural durabilities.

### **Functional Features:**

#### **1. Palpation of the prostate:**

- A. Normal prostate: Analog chestnuts size, transverse diameter 4cm, vertical diameter 3cm, anteroposterior diameter 2cm. In the middle there is a longitudinal shallow trench behind the prostate that is prostate ditch.
- B. Benign prostatic hyperplasia: prostatic degree I hyperplasia, prostate enlargement, simulated egg size, behind of prostate is flat, and the ditch is shallow in the middle.
- C. Benign prostatic hyperplasia: Prostate degree II hyperplasia, prostate enlargement, simulated the size of duck egg, the groove in the middle of prostate is disappeared.
- D. Benign prostatic hyperplasia: prostatic degree III hyperplasia, severe prostate enlargement, surface regularity, texture hard, simulated goose egg size, the bottom of the prostate can not touch.

#### **2. Rectal palpation:**

- A. Normal rectal
- B. Rectal polyps: Rectal posterior wall surface palpable nodules lumps, hard.

C. Early colorectal cancer: the posterior wall of the rectum surface palpable nodule mass, surface uneven, hard, it is development phase of colorectal cancer late stage.

### **3. Matters needing attention:**

A. should avoid prolonged exposure to sunlight to prevent deformation and discoloration.

B. when the model pollution, use gauze or cotton wool dipped in diluted cleaning wipe.

C. avoid models fall off, throw touch, any sharp objects scribing on the model or violence treat model will shorten the service life of the model.

**Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd.**

**Ada Med Supply Limited**

---

Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd. Building 3,  
No.1288 kungang Road, Xiaokunshan Town, Songjiang District, Shanghai, China

TEL: +86-13383897707

FAX: 021—57722933

Post code: 201612

E-mail: [edith@adahealthy.com](mailto:edith@adahealthy.com)

Website: [www.adahealthy.com](http://www.adahealthy.com)