

MANUAL

BIX-H70/3 Advanced Female Catheterization Model



Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd.

Ada Med Supply Limited

BIX-H70/3 Advanced Female Catheterization Model

1.Product Features:

The product advanced female catheterization model, is based on clinical and practice-based nursing care technical outline of the requirements ,newly developed by the company

The model is made of plastic material, made by die cast molding, injection, puncture, genitals, etc. are made of imported soft plastic.

With a vivid, real operation, easy accessibility, reasonable structure, durability and other characteristics.

2.Product Functions:

- Female Catheterization
- Female bladder lavage
- Ostomy drainage
- Intramuscular injection
- Femoral intramuscular
- Female Catheterization

Catheterization before the operation, intravesical water 200 ml firstly,

Hip abduction can be about 60 degrees, When the catheter is inserted into the bladder to about 4 cm, there is a real sense of obstruction.

Insert about 20 cm to the bladder. There are "urine" self discharge conduit.

三、 Usage

- 1.Medical University, School of Nursing, School of Health, clinical teaching in medical colleges and teaching students practical operation training
2. Hospital medical, nursing continuing education clinical teaching practice on training
- 3.Universal basic health units Clinical Training.

4.Maintenance Warranty

- 1.After an enema, catheterization operation training and other functions, to empty residue when operating
- 2.when not used for a long time, model should be wiped clean, packaged and placed large cool dry place, to extend the life.

Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd.

Ada Med Supply Limited

Shanghai Chinon Medical Model & Equipment Manufacturing Co.,Ltd. Building 3,
No.1288 kungang Road, Xiaokunshan Town, Songjiang District, Shanghai, China

Tel: +86-13383897707

Fax: 021—57722933

Post code: 201612

E-mail: edith@adahealthy.com

Website: www.adahealthy.com